




An ISO 13485-2016 & 2001 : 2015 certified company



Adult Diapers

- Non- woven one way top channels fluid into the absorbent pad & prevents it from flowing back
- Inner Absorbent pad
- Super Absorbent Polybeads turn fluid into gel. Users feel dry & comfortable
- Two pairs of refastenable tapes allow multiple adjustments and refastening on the frontal tape
- Multistrand Leg Elastic gives comfortable snug fit without pinching on skin
- Prevents leakage in the crotch area
- Unique maximum protection against leakage. Stops urine or faeces flowing towards the elastic
- Sizes: Small/ Medium/ Large/ Extra large
- Adult Diapers are an effective way to deal with the problem of elderly incontinence; the inability to exercise control over the bladder & the bowel
- They are a blessing for people suffering from Prostate or Urology disorders, Piles, Incontinence, Diabetes, or even patients who are bed-ridden & disabled

Underpads

- Soft Non- woven surface
- Ensures faster absorbency, Provides dry comfortable surface
- Unique design prevents leakage & ensure cleanliness
- Anti- bacterial super absorbency polymer ensures comfortable dry clean surface
- Available in Medium and Large sizes
- Incontinent patient's home care taking
- Post- maternity & after surgery use
- Hospital & nursing home care taking
- Also suitable for light discharges

Wrapping Sheets

- Effective barrier- Prevents microbial penetration and maintains sterility of surgical packs
- Penetrability- Allows steam and sterilizing gases to penetrate
- Aeration- Allows ethylene oxide to dissipate
- Minimal linting- During usage, linting is minimal
- Strength- Resists tears and punctures
- Drapeability- Conforms to a sterile pack and drapes will for ease of presentation
- Utility- Three weights and a range of sizes to meet your needs
- Ease-of-use- Easy to use by personnel

MFD & MKT by :

ramco

An ISO 13485-2016 & 2001 : 2015 certified company

RE-17/10 Raghunathpur Housing Estate
Raghunathpur, Kolkata - 700 059, WB, INDIA
Ph. - +91 9831 2322 47 / +91 90070 380 20

WORKSHOP
Panchanantala Industrial Area
Howrah, West Bengal, PIN-711 112