



Oxygen Mox Regulator

Features

- High flow rate against minimum pressure drops
- Accurate performance
- Simple design
- Easy installation
- Also available Double Stage Double Gauge Oxygen Regulator



Nitrous oxide Mox Regulator

Features

- High flow rate against minimum pressure drops
- Accurate performance
- Simple design
- Easy installation
- Also available Double Stage Double Gauge nitrous oxide Regulator



Double Stage Oxygen Regulator

Features

- Body: Forged brass with bright chrome plated
- Pressure gauge range 0 to 16 kg/cm², diameter- 2"
- Pressure gauge preset at 4.2 kg/cm²
- To protect over pressurization built in Safety valve provided
- Durable Neoprene rubber used
- To protect foreign matter Filter used at the center of inlet

Manufactured & Marketed by:

ramco

An ISO 13485:2003 & 9001:2008 certified company
RE-17/10 Raghunathpur Housing Estate
Raghunathpur, Kolkata - 700 059, WB, INDIA
Ph. - +91 9831 2322 47 / +91 90070 380 20

WORKSHOP

Panchanantala Industrial Area
Howrah, West Bengal, PIN-711 112